

news

uku – the flexible integrating sphere series

uku is the integrating spheres series of opsira.

uku series integrating spheres enable the fast and easy setup for measurements of the radiant power and the luminous flux of light sources as well as the reflection, the diffuse reflection and transmission of materials.

The idea is simple, exact and flexible. The uku series concept is based on planar surfaces allowing to adjust any type of detector, light source and material sample easily and without port adapters in a highly flexible way to the sphere ports.

A wide variety of accessories is available for the most flexible use possible. opsira offers various sizes of integrating spheres.

Please find further details here: www.opsira.com/uku



Pictures are ready for download here: www.opsira.com/downloads/press

Press contact: Ms Uta Vocke

vocke@opsira.de

www.opsira.de

January 14, 2014